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Zeng X, Z. S., Lu W, et al. (2020). Target Identification among Known Drugs by Deep Learning from Heterogeneous Networks. Chemical Science.

Jin J, X. N., Liu Y, et al. (2019). A novel S1P1 modulator IMMH002 ameliorates psoriasis in multiple animal models. Acta Pharmaceutica Sinica B.

Li H, S. J., Zhao M, et al. (2019). IL-17A can Promote Propionibacterium acnes-induced Sarcoidosis-like Granulomatosis in Mice. Frontiers in Immunology.

Lu L, Q. Y., Zhang X, et al. (2019). Shexiang Baoxin Pill Alleviates the Atherosclerotic Lesions in Mice via Improving Inflammation Response and Inhibiting Lipid Accumulation in the Arterial Wall. Mediators of Inflammation.

Zhou J, C. P., Li Z, et al. (2019). Gene delivery of TIPE2 attenuates collagen-induced arthritis by modulating inflammation. International Immunopharmacology.

Xu M, D. X. Y., Chen Q Y, et al. (2018). Effect of compound sophorae decoction on dextran sodium sulfate (DSS)-induced colitis in mice by regulating Th17/Treg cell balance. Biomedicine & Pharmacotherapy.

Pan Y, Z. Y., Wu H, et al. (2017). A Therapeutic Peptide Vaccine Against PCSK9. Scientific reports.

Hu X, D. Y., Chen X, et al. (2017). Immune Response of A Novel ATR-AP205-001 Conjugate Anti-hypertensive Vaccine. Scientific Reports.

Zhang A, W. K., Zhou C, et al. (2016). Knockout of microRNA-155 ameliorates the Th1/Th17 immune response and tissue injury in chronic rejection. The Journal of Heart and Lung Transplantation.

Yu M, L. Q., Li H H, et al. (2016). IL-9 Inhibits Viral Replication in Coxsackievirus B3-Induced Myocarditis. Frontiers in Immunology.

Wu Y, C. X., Ge X, et al. (2015). Isoliquiritigenin prevents the progression of psoriasis-like symptoms by inhibiting NF-κB and proinflammatory cytokines. Journal of Molecular Medicine.

Wei J, W. C., Wang M, et al. (2014). The GTPase-activating protein GIT2 protects against colitis by negatively regulating Toll-like receptor signaling. PNAS.

Junjie Wei, M. Z., Peixiao Liang, Yunfei Wei et al. (2013). Apolipoprotein E and its mimetic peptide suppress Th1 and Th17 responses in experimental autoimmune encephalomyelitis. Neurobiology of Disease.

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